8. Place enough amount of **BIOSEAL LIGHT CURED** along all pits and fissures extension, taking care to avoid the air entrapment;
9. Expose all faces of the sealed area to be light cured during 20 seconds.
10. Check occlusion and grind off any excess material, if necessary;
11. Check adhesion with an explorer.

CONTROL AND REAPPLICATION:

The applied sealant must be examined every 6 months. If it is totally absent, the reapplication must be done, following the use instructions. If parts of the sealant is present, use sharpened instruments to test its adherence. The sealant surface firmly adhered can receive a new sealant layer.

SPECIAL CARES:

Avoid direct contact of the product or of Etching Gel with the skin or mucous.

In case of contact, rinse thoroughly with water. DO NOT INGEST.

Protect eyes during polymerization procedure according to manufacturer instructions.

After using keep the syringes tightly closed to avoid exposure to the light.

Do not store **BIOSEAL LIGHT CURED** near other products containing Eugenol. Eugenol interferes in the polymerization.

Dispose the product in accordance with local regulations, ensuring its complete distortion, preventing its reuse and environment damages.

CONTRAINDICATION:

BIOSEAL LIGHT CURED should not be used in patients is known to have sensibility to any of the ingredients listed.


BIOSEAL LIGHT CURED is not recommended in pits and fissures with good congregation, in teeth with caries even initial and in teeth free of caries over 4 years.

WARNINGS:

Do not use material after expiry date.

The manufacture is not responsible for any damage resulting from a misuse or not provided for the instructions of use.
 "KEEP MATERIAL OUT OF THE REACH OF CHILDREN".

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or

1x5mL flask **BIOSEAL LIGHT CURED**
1x5mL flask Ataque Gel (37% Orthophosphoric Acid).
10 disposable brushes.

Reg. ANVISA: 10298550057.

COMPOSITION:

Bisphenol A Glycidylmethacrylate (35,6%); Methacrylate Groups; B.H.T.; Silicium Dioxide; Sodium Fluoride; Calcium Fluoride and Catalyst.

INDICATION:

BIOSEAL LIGHT CURED is indicated to seal pits and fissures, acting as a complement in caries prevention. It is indicated as a complement in caries prevention.

TECHNICAL INFORMATION:

BIOSEAL LIGHT CURED was developed to be applied on the occlusal surface of the teeth. After polymerization **BIOSEAL LIGHT CURED** forms a resistant and continuous film, used to seal pit and fissures preventing them against caries.

Using of sealants has been highly recommended considering its performance in reducing caries.

CHARACTERISTICS / PROPERTIES:

Light cured and radiopaque material; Excellent flow.

USE INSTRUCTIONS:

The efficiency of sealing depends on the correct penetration of the sealant in the conditioned enamel.

1. Clean occlusal surfaces including pits and fissures, with pumice powder and water;
2. Isolate the teeth to be sealed using cotton rolls or rubber dam and dry them well;
3. Etching with 37% Orthophosphoric Acid (a.e. Ataque Gel), during 30 seconds on the enamel;
4. Rinse with water during 30 seconds to remove the etching precipitates;
5. Dry thoroughly. The presence of humidity in the cleaned surface will impede the penetration of the resin in the conditioned area;
6. If the surface comes in contact with saliva, the etching and drying procedures must be repeated; (etch again for 10 seconds and dry);
7. The appropriately conditioned area should present a chalky white appearance; this not happening, repeat the procedure;